



Non Contact Forehead Thermometer (FR)

FR 1DZ1

■ Spec & Features



Non-Contact Forehead thermometer with Silent Glow™

- Probe no need to contact forehead
- Reference mode selection (rectal, oral, axillary)
- Offset for personalized readings
- Object Mode & Body Mode switch able
- 30 memories in memory mode
- °C / °F Switchable
- Silent glow
- Measuring range : 34 °C~42.2 °C Body mode
0~100 °C(Object mode)
- Low battery indication
- Auto power-off / 60"